




- ✓ Selective Blue Light Blocking (blue-violet: 380 to 450 nm)
- ✓ Class 1 UV Filter (UVA: 316 to 380 nm and UVB: 280 to 315 nm)
- ✓ Available in Spheric, Toric, Multifocal and Multifocal Toric geometries
- ✓ 5.3 billion parameter combinations to match any patient's eye and vision needs

CHARACTERISTICS	MATERIAL	GEOMETRIES	ADVANCED OPTICAL DESIGNS
<ul style="list-style-type: none"> ■ Class 1 UV Filter ■ Blue Light Filter ■ Monthly Replacement ■ Custom Lens 	<ul style="list-style-type: none"> ■ Silicone Hydrogel ■ 75% H₂O ■ Dk 60 ■ Low Dehydration ■ Handling Tint, Green 	<ul style="list-style-type: none"> ■ Spheric ■ Toric ■ Multifocal ■ Multifocal toric 	Multifocal - Design CD - CN

BLU:GEN & BLU:KIDZ FITTING PROTOCOL

Ø (mm)	11.50	12.00	12.50	13.00	13.50	14.00	14.50	15.00	15.50	16.00	16.50
Base Curve (mm)	6.50 - 8.30	6.50 - 8.30	6.50 - 8.60	6.50 - 8.90	6.80 - 9.20	7.10 - 9.50	7.40 - 9.80	7.70 - 9.80	8.00 - 9.80	8.30 - 9.80	8.60 - 9.80
Fitting Rule Km = (K1 + K2)/2	Km +0.0			Km +0.1	Km +0.3	Km +0.5	Km +0.7	Km +0.9	Km +1.1	Km +1.3	

PARAMETER RANGE




■ Standard lens ■ Blu:gen/Blu:kidz


Inspire your eyes

markennovy

THE markennovy DIFFERENCE



An unrivalled combination of lens parameters, geometries and materials



Driving loyalty through enhancing the patient experience leading to better practice profitability



Fast delivery – diagnostic lenses shipped in 72hours



Technical support staffed by opticians, for opticians



Solely focused on eye care professionals, protecting your business against the internet

Contact your **Business Development Manager** for details

Customer Care: T. 0800 328 0610 F. 0800 328 0649
International Customer Care: T. +34 916 496 123
mkservices@markennovy.com