

WEEKLY REPLACEMENT

SEVEN RX

HYDROGEL



Seven RX is our weekly lens developed from our proven hydrogel material. This contact lens is lathed to offer a wider range of toric parameters than its moulded version.







PARAMETERS

BASE CURVE (mm)	SPH, MF : 8.50 TOR, MFT : 8.70
DIAMETER (mm)	14.40
SPHERES (D)	SPH ± 30.00 (0.50 after ± 6.00) TOR, MF, MFT ± 30.00 (0.50 after +4.00/-6.00)
CYLINDERS (D)	-0.75 to -7.75 (0.50)
AXES (°)	All (5°)
ADDITION	CN +1.25/+2.25 CD +1.50/+2.50


MATERIAL

TYPE	Filcon 4 (19) [55%]
WATER CONTENT	55%
CENTRAL THICKNESS (-3,00 dpt)	0.10
DK (ISO 9913-1-1998)	19
DK/T (-3.00 D)	19
HANDLING TINT	BLUE
UV BLOCKING	NO
PACK SIZE	12 Lenses
MANUFACTURING PROCESS	LATHED

OPTICAL DESIGN



Lens design parameters may change depending on the power



Inspire your eyes
markennovy

**CAPRICORNIA
CONTACT LENS**

CONTACT CHRIS HAROUS:
T: 0400 217 007
e-mail: chris.harous@markennovy.com

Telephone: +61 (7) 3208 8500
Tollfree (Australia): 1300 650 994
Fax: +61 (7) 3208 8506
info@capcl.com.au

markennovy