


MAANDLENS

XTENSA RX

HYDROGEL







PARAMETERS

BASIS CURVE (mm)	SFE, MF : 8.50 TOR, MFT : 8.70
DIAMETER (mm)	14.40
SFERISCHE STERKTES (D)	SFE ± 30.00 (0.50 na ± 6.00) TOR, MF, MFT ± 30.00 (0.50 na $+4.00/-6.00$)
CILINDER (D)	-0.75 to -7.75 (0.50)
AS ($^{\circ}$)	Alle (5°)
ADDITIE	CN +1.25/+2.25 CD +1.50/+2.50

MATERIAAL

TYPE	Filcon 4 (19) [55%]
WATERGEHALTE	55%
CENTRALE DIKTE (-3.00D)	0.10 (SFE)
DK (ISO 9913-1-1998)	19
DK/T (-3.00D)	22
HANDLINGS TINT	BLAUW
UV FILTER	NEE
VERPAKKINGS GROOTTE	12 Pack
PRODUCTIEPROCES	LATHED

LENSONTWERPEN



De parameters van het lensontwerp kunnen afhankelijk van de sterkte variëren.

NEEM VOOR MEER INFORMATIE CONTACT OP MET UW ACCOUNT MANAGER OF BEL:

Customer Care NL: 0800 023 2029
Customer Care BE: 0800 799 18
E: klantenservice@markennovy.com

markennovy