



Versa One Day Toric is a made-to-order daily disposable contact lens. Available in **65 spheres, 15 cylinders and 36 axes**, Versa gives you a wider range of power and axis options to offer your patients unmatched astigmatic correction in a daily replacement. Its proven lens design and hydrogel material provide patients with stable vision and excellent comfort all day, every day.

REPLACEMENT	Daily
PACK SIZE	30 & 90 lenses
GEOMETRY	🕒 Toric
MATERIAL CHARACTERISTICS	Hydrogel (55% H ₂ O)
	Dk 19
	MPa 0.47
	Handling tint Yes
PARAMETERS	BC (mm) 8.70
	Ø (mm) 14.40
	SPH (D) ±10.00 (0.50D steps after ±6.00)
	CYL (D) -0.75 to -5.75 (0.50D steps after -2.75)
	AXIS (°) 0 to 180 in 5° steps


CONTACT YOUR BUSINESS DEVELOPMENT MANAGER FOR MORE DETAILS

Customer Care: T. 0800 328 0610 F. 0800 328 0649
 International Customer Care: T. +34 916 496 123
 mkservices@markennovy.com

OUR DIFFERENCE



Exclusive focus on the Eye Care Professional




Exceptional combination of parameters, geometries and materials



Optician-led Customer Care



Individually Yours



Fast delivery – diagnostic lenses shipped in 72 hours

CONTACT YOUR BUSINESS DEVELOPMENT MANAGER FOR MORE DETAILS

Customer Care: T. 0800 328 0610 F. 0800 328 0649
International Customer Care: T. +34 916 496 123
mkservices@markennovy.com

Inspire your eyes
markennovy