

MONTHLY REPLACEMENT

XTENSA

HYDROGEL




Xtenso is a hydrogel lens developed from our proven hydrogel material. It offers a fundamental range of parameters to meet standard spherical, toric, and multifocal prescriptions. Its toric range is complemented by Xtenso Rx for patients whose astigmatic needs extend beyond standard.



PARAMETERS	
BASE CURVES (mm)	SPH, TOR 8.70 MF 8.60
DIAMETERS (mm)	SPH, TOR 14.40 MF 14.30
SPHERES (D)	SPH -10.00 to +8.00 (0.50 after ± 6.00) TOR -8.00 to +6.00 (0.50 after +4.00/-6.00) MF ± 10.00 (0.50 after +4.00/-6.00)
CYLINDERS (D)	-0.75 to -2.75 (0.50)
AXES (°)	0 to 180 (10°)
ADDITION	CN +1.25/+2.25 CD +1.50/+2.50

MATERIAL	
TYPE	Filcon 4 (19) [55%]
WATER CONTENT	55%
CENTRAL THICKNESS (- 3,00 dpt)	0.09 (ASP) 0.10 (TOR, CYL -1.75) 0.13 (MF-CD) 0.12 (MF-CN) 0.11 (MFT-CD) 0.14 (MFT-CN)
DK (ISO 9913-1-1998)	19
DK/T (-3.00 D)	21 (ASP) 19 (TOR) 15 (MF-CD) 16 (MF-CN) 17 (MFT-CD) 14 (MFT-CN)
HANDLING TINT	BLUE
UV BLOCKING	NO
PACK SIZE	6 Lenses
MANUFACTURING PROCESS	MOULDED



Lens design parameters may change depending on the power

CONTACT YOUR BUSINESS DEVELOPMENT MANAGER FOR DETAILS

Customer Care: T. 0800 328 0610 F. 0800 328 0649
International Customer Care: T. +34 916 496 123
mkservices@markennovy.com

markennovy