

# 3-MONTHLY REPLACEMENT

## EQUILIBRIA

HYDROGEL



**Equilibria** provides a non-silicone option, featuring good water retention and tensile properties, for patients already accustomed to a 3-monthly lens replacement.

-  SPHERIC
-  MULTIFOCAL
-  TORIC
-  MULTIFOCAL TORIC

### PARAMETERS

<b>BASE CURVES (mm)</b>	7.70 to 9.80 (0.30)
<b>DIAMETER (mm)</b>	14.50
<b>SPHERES (D)</b>	SPH, TOR: $\pm 30.00$ (0.25)
	MF, MFT: $\pm 23.00$ (0.25)
<b>CYLINDERS (D)</b>	-0.75 to -8.00 (0.25)
<b>AXES (°)</b>	All (5°)
<b>ADDITIONS</b>	1.00 to 3.00 (0.50) CD/CN


### MATERIALE

<b>TYPE</b>	Filcon 2 (24) [59%]
<b>WATER CONTENT</b>	59%
<b>DK (ISO 9913-1-1998)</b>	24
<b>HANDLING TINT</b>	BLUE
<b>MODULUS</b>	0.31
<b>PACK SIZE</b>	Single and 2-pack
<b>MANUFACTURING PROCESS</b>	LATHED
<b>CoF</b>	1.17

### Calculate your lens

<b>LØ (mm)</b>	14.50
<b>CB (mm)</b>	7.70 - 9.80
<b>FITTING RULE</b>	0.8
<b><math>K_m = (K_1 + K_2) / 2</math></b>	

### OPTICAL DESIGN



Lens design parameters may change depending on the power

### CONTACT YOUR BUSINESS DEVELOPMENT MANAGER FOR DETAILS

Customer Care: T. 0800 328 0610 F. 0800 328 0649  
 International Customer Care: T. +34 916 496 123  
 mkservices@markennovy.com

**markennovy**