

MONTHLY REPLACEMENT

XTENSA RX





HYDROGEL



PARAMETERS	
BASE CURVE (mm)	SPH, MF : 8.50 TOR, MFT : 8.70
DIAMETER (mm)	14.40
SPHERES (D)	SPH ± 30.00 (0.50 after ± 6.00) TOR, MF, MFT ± 30.00 (0.50 after +4.00/-6.00)
CYLINDERS (D)	-0.75 to -7.75 (0.50)
AXES (°)	All (5°)
ADDITION	CN +1.25/+2.25 CD +1.50/+2.50

MATERIAL	
TYPE	Filcon 4 (19) [55%]
WATER CONTENT	55%
CENTRAL THICKNESS (-3.00 D)	0.10 (SPH)
DK (ISO 9913-1-1998)	19
DK/T (-3.00 D)	19 (SPH)
HANDLING TINT	BLUE
UV BLOCKING	NO
PACK SIZE	6 Lenses
MANUFACTURING PROCESS	LATHED

OPTICAL DESIGN



Lens design parameters may change depending on the power

CONTACT YOUR BUSINESS DEVELOPMENT MANAGER FOR DETAILS

Customer Care: T. 0800 328 0610 F. 0800 328 0649
International Customer Care: T. +34 916 496 123
mkservices@markennovy.com

markennovy